


OSSIOfiber[®] Cannulated Trimmable Fixation Nail System Surgical Technique

The OSSIOfiber[®] Cannulated Trimmable Fixation Nail System includes strong and biointegrative trimmable nails conveniently packaged with corresponding sterile, disposable instrumentation. The system provides a solution to securely fixate bone fractures or osteotomies for natural bone healing, ultimately leaving the patient renewed without permanent metal hardware.



For 3.0 x 50mm: two 1.1mm K-wire, 3.0 cannulated drill bit, depth gauge and tamp For 4.0x70mm: three 1.4mm K-wire, 4.0 cannulated drill bit, depth gauge and tamp

DOC-0001765 Rev3 1/2022


® OSSIO and OSSIOfiber are registered trademarks of OSSIO Ltd. All rights reserved. OSSIO Inc. 300 Tradecenter Drive, Suite 3690, Woburn, MA 01801 For more on OSSIO and OSSIOfiber[®], please visit ossio.io or call 833-781-7373

Not available for sale outside of the US, Speak to your local sales representative for product availability.

Sterile, disposable instrumentation

Convenient and streamlined OR preparation

Surgical Technique



Step 1: Place K-Wires

Insert K-Wires in the desired location.



Step 4: Tamp nail into tunnel


Place tapered end of nail over k-wire followed by the cannulated tamp. Use a mallet to tap the proximal side of the tamp to advance the implant into the tunnel.



Step 2: Measure depth

Place provided depth gauge over k-wire to measure the tunnel length. Mark desired length on nail with a surgical marker.

Note: Pre-trimming with an oscillating saw is an option at this step.



Step 3: Prepare tunnel


Place the cannulated drill over the k-wire to create pilot hole to measured nail depth.

5



Step 5: Remove instrumentation

Remove K-wire to leave the implant after final fixation is achieved. Trim the nail flush to bone with an oscillating saw (if not pre-trimmed in step 2).



Note: If recessed placement is desired, turn the tamp knob to expose the inner tip of the tamp. Use a mallet to drive the nail into the hole just below the bone surface.