



OSSIOfiber™

Hammertoe Fixation System


The Strong and Bio-Integrative Implant

Comprised of **OSSIOfiber™ Intelligent Bone Regeneration Technology**

- Strong, secure fixation for PIP fusions
- Fully integrates into native anatomy without adverse inflammation, as proven in pre-clinical studies
- Visible yet artifact-free on CT and X-ray, MRI safe

Three sizes
(2.5 mm, 2.9 mm, and 3.2 mm)

Two orientations
(0° and 10°)



Convenient Disposable Kit

- Disposable K-Wire, Drill Bit, and Inserter provided
- Implants packaged with disposable instruments and as stand-alone for cost savings
- No additional handling requirements or freezer storage required





Kit ordering information

Disposable kits (implant + instrumentation) and stand-alone implants are available in 3 sizes (small, medium, and large) and in 2 orientations (straight and angled).

| Size | Description | Product Code |
|--------|---|--------------|
| Small | OSSIOfiber™ Hammertoe Fixation Implant 2.5 x 16 mm, 0° Straight* | OF10025160 |
| | OSSIOfiber™ Hammertoe Fixation System 2.5 x 16 mm, 0° Straight Implant and Instrument Kit† | OF20025160 |
| | OSSIOfiber™ Hammertoe Fixation Implant 2.5 x 16 mm, 10° Angled* | OF10025161 |
| | OSSIOfiber™ Hammertoe Fixation System 2.5 x 16 mm, 10° Angled Implant and Instrument Kit† | OF20025161 |
| Medium | OSSIOfiber™ Hammertoe Fixation Implant 2.9 x 19 mm, 0° Straight* | OF10029190 |
| | OSSIOfiber™ Hammertoe Fixation System 2.9 x 19 mm, 0° Straight Implant and Instrument Kit† | OF20029190 |
| | OSSIOfiber™ Hammertoe Fixation Implant 2.9 x 19 mm, 10° Angled* | OF10029191 |
| | OSSIOfiber™ Hammertoe Fixation System 2.9 x 19 mm, 10° Angled Implant and Instrument Kit† | OF20029191 |
| Large | OSSIOfiber™ Hammertoe Fixation Implant 3.2 x 21 mm, 0° Straight* | OF10032210 |
| | OSSIOfiber™ Hammertoe Fixation System 3.2 x 21 mm, 0° Straight Implant and Instrument Kit† | OF20032210 |
| | OSSIOfiber™ Hammertoe Fixation Implant 3.2 x 21 mm, 10° Angled* | OF10032211 |
| | OSSIOfiber™ Hammertoe Fixation System 3.2 x 21 mm, 10° Angled Implant and Instrument Kit† | OF20032211 |

*Implants are packaged pre-loaded on an inserter.

†Instrument kits include a Drill Bit and K-Wire.

