

OSSIOfiber[®]

Threaded Trimmable Fixation Nail

Intentionally designed for
Metacarpal Fracture Fixation

Measure



Trim



Transform




Thousands of **continuous natural mineral fibers** are organized into a custom matrix designed to optimize **strength** and enable safe **bio-integration**.

Unmatched Rotational Stability

Outperforming all major metal competitors, the OSSIOfiber[®] Threaded Trimmable Fixation Nail provides greater rotational stability compared to widely-used metal competitors¹.

One Size Fits All

The nail length can be customized intraoperatively, reducing the need for extensive inventory.



Avoid Nightmare Revision Surgeries

OSSIOfiber[®] bio-integrative technology allows for drilling through the existing nail, obviating the need for difficult removal of intramedullary nails and screws from the metacarpal.

Hand & Wrist Threaded Trimmable Fixation Nail Ordering Info

The OSSIOfiber® Threaded Trimmable Fixation Nail System includes one instrumentation kit and one implant.



| Description | Product Code |
|---|--------------|
| OSSIOfiber® Threaded Trimmable Fixation Nail, 4.5 X 70 mm | OF1084570S |
| OSSIOfiber® Threaded Trimmable Fixation Nail, 4.5 X 70 mm System, H/W | OF2084510S |
| OSSIOfiber® Threaded Trimmable Fixation Nail Sizing Guide Pack | THN00480 |

Hand & Wrist Sterilized and Disposable Instrument Pack

- 1.1 mm K-wire (x2)
- Cannulated Drill Bit
- Implant Sharpener
- Adjustable Driver



Sizing Guides

- 4.5 mm 
- 5.5 mm 

Indications

The OSSIOfiber® Threaded Trimmable Fixation Nails are indicated for maintenance of alignment and fixation of bone fractures, comminuted fractures, fragments, osteotomies, arthrodesis, and bone grafts, of the upper extremity, fibula, knee, ankle and foot in the presence of appropriate brace and/or immobilization **in adults and children (2-12 years) and adolescents (12-21 years)** in which growth plates have fused or in which growth plates will not be crossed by fixation.

1. Data on file at OSSIO

DOC-0003664 09/2024

OSSIO® and OSSIOfiber® are registered trademarks of OSSIO Ltd. All rights reserved 2024.

Refer to the product Instructions for Use for warnings, precautions, indications, contraindications and technique.

Some products mentioned in this document may not be currently available or approved for sale in your country.

Speak to your local sales representative/distributor for product availability.

300 Tradecenter Drive, Suite 3690, Woburn, MA 01801

Visit [OSSIO.com](https://www.ossio.com)



OSSIO[®]
Restore. Regrow. Renew.